

GOAT ANTI-SMAD6 / MADH6 ANTIBODY

SKU: EB05150



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	MAD (mothers against decapentaplegic, Drosophila) homolog 6 Mothers against decapentaplegic, drosophila, homolog of, 6 SMAD6 MADH7 MAD, mothers against decapentaplegic homolog 6 (Drosophila) MADH6
Blocking Peptide	EBP05150
Immunogen	Peptide with sequence PCWLEILLNNPR, from the C Terminus of the protein sequence according to .
Peptide Sequence	PCWLEILLNNPR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	4091
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 55kDa band seen in Jurkat lysate. Recommended for use at 1-3µg/ml.