

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SMURF1 ANTIBODY

SKU: EB05268

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias E3 ubiquitin ligase SMURF1|SMURF1

Names

Accession AAF08298 ID

Blocking

EBP05268 Peptide

Immunogen AAF08298. Peptide with sequence DVRGSLQTPQNRPHC, from the internal region of the protein sequence according to

Peptide

DVRGSLQTPQNRPHC Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

Gene ID

57154

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Western Western blot: No signal obtained yet but low background observed in HeLa and NIH-3T3 cell lysates at upto

Blot $2\mu g/ml$.