

GOAT ANTI-SMURF1 ANTIBODY

SKU: EB05268

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	E3 ubiquitin ligase SMURF1 SMURF1
Accession ID	AAF08298
Blocking Peptide	EBP05268
Immunogen	Peptide with sequence DVRGSLQTPQNRPHC, from the internal region of the protein sequence according to AAF08298.
Peptide Sequence	DVRGSLQTPQNRPHC
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	57154
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in HeLa and NIH-3T3 cell lysates at upto 2µg/ml.