

GOAT ANTI-SMURF1 ANTIBODY

SKU: EB06975



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	SMURF1 SMAD specific E3 ubiquitin protein ligase 1 HGNC:16807 KIAA1625 E3 ubiquitin ligase SMURF1 Smad ubiquitination regulatory factor 1 Smad-specific E3 ubiquitin ligase 1
Accession ID	NP_065162.1; NP_851994.1; NP_001186776.1
Blocking Peptide	EBP06975
Immunogen	Peptide with sequence C-PSEGSLEDEELPAQ, from the internal region of the protein sequence according to NP_065162.1; NP_851994.1; NP_001186776.1.
Product Comments	This antibody is expected to recognise all 3 reported isoforms
Peptide Sequence	C-PSEGSLEDEELPAQ
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	57154
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

