

GOAT ANTI-SMURF2 ANTIBODY

SKU: EB05272



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	E3 ubiquitin ligase SMURF2 SMURF2
Accession ID	NP_073576
Immunogen	Peptide with sequence SNPGRRNGPVKLR-C, from the N Terminus of the protein sequence according to NP_073576.
Peptide Sequence	SNPGRRNGPVKLR-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	64750
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: Three consistent bands observed in Human Brain extracts at approx 80, 50 and 40kDa with 1 hour incubation of primary.