

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SMURF2 ANTIBODY

SKU: EB05272

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias E3 ubiquitin ligase SMURF2|SMURF2

Names

Accession NP 073576

Peptide with sequence SNPGGRRNGPVKLR-C, from the N Terminus of the protein sequence according to

Immunogen NP_073576.

Peptide

SNPGGRRNGPVKLR-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse **Species**

Human

64750 **Gene ID**

Product Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/aspiring_medium.png$

Western Western blot: Three consistent bands observed in Human Brain extracts at approx 80, 50 and 40kDa with 1 hour

Blot incubation of primary.