

GOAT ANTI-SNURPORTIN1 ANTIBODY

SKU: EB06245



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	SNUPN snurportin1 snuportin-1 RNA, U transporter 1 SNURPORTIN-1 KPNBL RNUT1
Accession ID	NP_005692
Immunogen	Peptide with sequence SHSPDHPGCLMEN, from the C Terminus of the protein sequence according to NP_005692.
Peptide Sequence	SHSPDHPGCLMEN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	10073
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Western Blot: Preliminary experiments gave no signal but low background in Hela and Human Brain extracts at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA