

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SNURPORTIN1 ANTIBODY

SKU: EB06245

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias SNUPN|snurportin1|snuportin-1|RNA, U transporter 1|SNURPORTIN-1|KPNBL|RNUT1

Names

Accession NP 005692

Immunogen Peptide with sequence SHSPDHPGCLMEN, from the C Terminus of the protein sequence according to NP_005692.

Peptide SHSPDHPGCLMEN

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted Human

Species

Human 10073

Gene ID

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:2000.

Limit

Blot

Western Blot: Preliminary experiments gave no signal but low background in Hela and Human Brain extracts at up Western to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?

Application

Pep-ELISA **Type**