

GOAT ANTI-SNX10 ANTIBODY

SKU: EB05891



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	MGC33054 sorting nexin 10 SNX10
Accession ID	NP_037454.2
Blocking Peptide	EBP05891
Immunogen	Peptide with sequence SHGCKVNTAPQES, from the C Terminus of the protein sequence according to NP_037454.2.
Peptide Sequence	SHGCKVNTAPQES
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	29887
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Western Blot: No signal obtained yet but low background observed in Human Liver, Human Testis, 293 and Raji lysates at upto 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA