

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-SNX10 ANTIBODY**

**SKU:** EB05891

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias MGC33054|sorting nexin 10|SNX10

**Names** 

Accession NP 037454.2 ID

**Blocking** 

Sequence

EBP05891 **Peptide** 

Immunogen Peptide with sequence SHGCKVNTAPQES, from the C Terminus of the protein sequence according to NP 037454.2.

**Peptide SHGCKVNTAPQES** 

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** Refrigerated

Instructions

**Predicted** 

Human **Species** 

Human Gene ID

29887

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:16000.

Limit

Western Blot: No signal obtained yet but low background observed in Human Liver, Human Testis, 293 and Raji Western Blot

lysates at upto 1µg/ml. We would appreciate any feedback from people in the field - have any results been

reported with other antibodies/lysates?

Application

Pep-ELISA **Type**