

GOAT ANTI-SNX11 ANTIBODY

SKU: EB05892



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	MGC111019 sorting nexin 11 SNX11
Names	
Accession ID	NP_037455.2; NP_689450.1
Blocking Peptide	EBP05892
Immunogen	Peptide with sequence C-PLDPGQLETVLEK, from the C Terminus of the protein sequence according to NP_037455.2; NP_689450.1.
Product Comments	Both variants (NP_037455.2; NP_689450.1) encode the same protein.
Peptide Sequence	C-PLDPGQLETVLEK
Shipping Instructions	Refrigerated
Predicted Species	Human, Cow
Human Gene ID	29916
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	No signal obtained yet but low background observed in Human Kidney, Human Testis, 293 and HepG2 lysates at upto 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA