


Telephone: (650) 697-3600

GOAT ANTI-SNX8 ANTIBODY

SKU: EB06219



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

SNX8|sorting nexin 8|Mvp1 **Alias**

Names

Accession NP_037453.1 ID

Blocking

EBP06219 **Peptide**

Immunogen

Peptide with sequence TPPCSPPEDGLCPH, from the C Terminus of the protein sequence according to

NP_037453.1.

Peptide

TPPCSPPEDGLCPH

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted

Human

Species Reactive

Species

Human

Human

29886

Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results

In paraffin embedded Human Kidney shows vesicular staining in the cytoplasm of preferentially DCT epithelial

cells. Recommended concentration, 4-8µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:8000.

Limit

Type

Application Pep-ELISA, IHC

GALLERY IMAGES

