


Telephone: (650) 697-3600

GOAT ANTI-SORL1 / LR11 (C TERMINUS) ANTIBODY

SKU: EB05045



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage Instructions

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

sortilin-related receptor containing LDLR class A repeats|FLJ39258|FLJ21930|C11orf32|mosaic protein LR11|SorLA-1|gp250|SORLA|LRP9|sortilin-related receptor, L(DLR class) A repeats-containing|LR11|SORL1

Names Accession

Alias

ID

NP_003096.1

Blocking

EBP05045

Peptide Immunogen

Peptide with sequence C-TGFSDDVPMVIA, from the C Terminus of the protein sequence according to

NP 003096.1.

Peptide

C-TGFSDDVPMVIA

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method

using the immunizing peptide.





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Shipping Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog, Cow

Species Reactive

Human **Species**

Human

Gene ID

6653

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results Paraffin embedded Human Brain (Cerebellum and Cortex). Recommended concentration: 5-7µg/ml.

ELISA

Detection Antibody detection limit dilution 1:64000.

Limit

Type


Application

Pep-ELISA, IHC


SELECTED REFERENCES

[{"pmid": 21040701, "intro": "This antibody has been successfully used in Mouse in the following paper:", "title": "Implication of mouse Vps26b-Vps29-Vps35 retromer complex in sortilin trafficking.", "author": "Kim E, Lee Y, Lee HJ, Kim JS, Song BS, Huh JW, Lee SR, Kim SU, Kim SH, Hong Y, Shim I, Chang KT.", "journal": "Biochem Biophys Res Commun. 2010 Dec 10;403(2):167-71."}]

GALLERY IMAGES











Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

