

GOAT ANTI-SOX11 ANTIBODY

SKU: EB11617



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	SRY (sex determining region Y)-box 11 Sox11
Accession ID	NP_003099.1
Blocking Peptide	EBP11617
Immunogen	Peptide with sequence C-DWLEANFSDLVFTY, from the C Terminus of the protein sequence according to NP_003099.1.
Peptide Sequence	C-DWLEANFSDLVFTY
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Pig, Cow, Xenopus, Danio rerio
Human Gene ID	6664
Mouse Gene ID	20666
Rat Gene ID	84046
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)