Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-SOX11 ANTIBODY**

**SKU:** EB11617

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias SRY (sex determining region Y)-box 11|Sox11

**Names** 

Accession NP 003099.1 ID

**Blocking** 

EBP11617 **Peptide** 

Peptide with sequence C-DWLEANFSDLVFTY, from the C Terminus of the protein sequence according to Immunogen

NP 003099.1.

**Peptide** C-DWLEANFSDLVFTY Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat, Pig, Cow, Xenopus, Danio rerio **Species** 

Human

6664 Gene ID

Mouse **Gene ID** 

20666

Rat Gene ID 84046

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

**ELISA** 

**Detection** 

Antibody detection limit dilution 1:8000.

Limit

**Application** Pep-ELISA

**Type**