

GOAT ANTI-SOX11 / SRY-BOX11 ANTIBODY

SKU: EB07386



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	transcription factor SOX-11 SRY-box 11 SRY-related HMG-box gene 11 SRY (sex-determining region Y)-box 11 SRY (sex determining region Y)-box 11 SOX11
Accession ID	NP_003099.1
Blocking Peptide	EBP07386
Immunogen	Peptide with sequence KTSKGSSKKCGK, from the internal region of the protein sequence according to NP_003099.1.
Peptide Sequence	KTSKGSSKKCGK
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	6664
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Preliminary experiments in Human Testis lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA