

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SOX8 (AA 203 TO 217) ANTIBODY

SKU: EB08216

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias MGC24837|SRY (sex determining region Y)-box 8

Names

Accession NP 055402.2 ID

Blocking

EBP08216 **Peptide**

Immunogen NP_055402.2. Peptide with sequence C-KAEAGLGDGHHHGDH, from the internal region of the protein sequence according to

Peptide C-KAEAGLGDGHHHGDH Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

30812 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:16000.

Limit

Application Pep-ELISA **Type**