

GOAT ANTI-SOX8 (AA 203 TO 217) ANTIBODY

SKU: EB08216



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	MGC24837 SRY (sex determining region Y)-box 8
Accession ID	NP_055402.2
Blocking Peptide	EBP08216
Immunogen	Peptide with sequence C-KAEAGLGDGHHHGDH, from the internal region of the protein sequence according to NP_055402.2.
Peptide Sequence	C-KAEAGLGDGHHHGDH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	30812
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)