

Telephone: (650) 697-3600

GOAT ANTI-SOX8 (AA 247 TO 258) ANTIBODY

SKU: EB08185

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

e erest

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias SOX8|MGC24837|SRY (sex determining region Y)-box 8

Names

Accession NP 055402.2 ID

Blocking

EBP08185 **Peptide**

Peptide with sequence C-EGRRPVDSGRQN, from the internal region of the protein sequence according to Immunogen

NP 055402.2.

Peptide C-EGRRPVDSGRQN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat **Species**

Human

30812 Gene ID

Mouse

20681

Gene ID

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:8000.

Limit

Application Pep-ELISA **Type**

> Goat Anti-SOX8 (aa 247 to 258) Antibody https://everestbiotech.com/products/goat-anti-sox8-aa-247-to-258-antibody/