

GOAT ANTI-SOX8 (AA 247 TO 258) ANTIBODY

SKU: EB08185



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. |
| Unit Size | 100 µg |
| Storage Instructions | Aliquot and store at -20°C. Minimize freezing and thawing. |
| Synonym / Alias Names | SOX8 MGC24837 SRY (sex determining region Y)-box 8 |
| Accession ID | NP_055402.2 |
| Blocking Peptide | EBP08185 |
| Immunogen | Peptide with sequence C-EGRRPVDSGRQN, from the internal region of the protein sequence according to NP_055402.2. |
| Peptide Sequence | C-EGRRPVDSGRQN |
| Purification Method | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| Shipping Instructions | Refrigerated |
| Predicted Species | Human, Mouse, Rat |
| Human Gene ID | 30812 |
| Mouse Gene ID | 20681 |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png |
| ELISA Detection Limit | Antibody detection limit dilution 1:8000. |
| Application Type | Pep-ELISA |

DOCUMENTS

- [Data Sheet](#)