


GOAT ANTI-SOX9B (ZEBRAFISH AA5-18) ANTIBODY

SKU: EB11444



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias sox9b| wu:fb18d01| SRY-box containing gene 9b| fb18d01

Names

Accession ID NP_571719.1

Blocking Peptide EBP11444

Immunogen Peptide with sequence QRGLKMSVSGAPSP-C, from the N Terminus of the protein sequence according to NP_571719.1.

Product Comments No cross-reactivity with Sox9a

Peptide Sequence QRGLKMSVSGAPSP-C

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Zebrafish
Reactive Species	Human
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Prostate. Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

