

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: <u>(650)</u> 697-3600

## **GOAT ANTI-SP1(MOUSE) ANTIBODY**

**SKU:** EB09068



Telephone: (650) 697-3600



250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Sp1-1|AI845540|AA450830|1110003E12Rik|trans-acting transcription factor 1|Sp1

Names

**Accession** NP\_038700.2

**Blocking** 

EBP09068 **Peptide** 

Peptide with sequence C-DGIHDETAGGEGEN, from the internal region of the protein sequence according to **Immunogen** 

NP 038700.2.

**Peptide** Sequence

C-DGIHDETAGGEGEN

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Mouse, Rat, Pig **Species** 

Reactive

Mouse **Species** 

Mouse 20683

Gene ID

**Product** 

Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:8000.

Limit

Western

Approx 100kDa band observed in nuclear lysates of cell line NIH3T3 (calculated MW of 80.4kDa according to NP\_038700.2). This molecular weight is routinely observed by other sources. Recommended concentration:

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**Blot**  $0.3-1\mu g/ml$ .

**Application** 

Pep-ELISA, WB

**Type** 

## **GALLERY IMAGES**



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa