

GOAT ANTI-SP1(MOUSE) ANTIBODY

SKU: EB09068



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Sp1-1 AI845540 AA450830 1110003E12Rik trans-acting transcription factor 1 Sp1
Accession ID	NP_038700.2
Blocking Peptide	EBP09068
Immunogen	Peptide with sequence C-DGIHDETAGGEGEN, from the internal region of the protein sequence according to NP_038700.2.
Peptide Sequence	C-DGIHDETAGGEGEN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Mouse, Rat, Pig
Reactive Species	Mouse
Mouse Gene ID	20683
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 100kDa band observed in nuclear lysates of cell line NIH3T3 (calculated MW of 80.4kDa according to NP_038700.2). This molecular weight is routinely observed by other sources. Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

