


GOAT ANTI-SPAG16 ANTIBODY

SKU: EB09825



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names SPAG16|WDR29|WD repeat domain 29|sperm-associated WD repeat protein|sperm associated antigen 16|PF20|MGC87036|FLJ37717|FLJ22724|DKFZp666P1710

Accession ID NP_078808.3

Blocking Peptide EBP09825

Immunogen Peptide with sequence C-KLHYASYEPTIR, from the internal region of the protein sequence according to NP_078808.3.

Product Comments This antibody is expected to recognize isoform 1 (NP_078808.3) only.

Peptide Sequence C-KLHYASYEPTIR

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Dog, Cow
Reactive Species	Human
Human Gene ID	79582
Mouse Gene ID	66722
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Testis. Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

