


Telephone: (650) 697-3600

GOAT ANTI-SPAG16 ANTIBODY

SKU: EB09825



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym / **Alias**

SPAG16|WDR29|WD repeat domain 29|sperm-associated WD repeat protein|sperm associated antigen 16|PF20|MGC87036|FLJ37717|FLJ22724|DKFZp666P1710

Names

Accession

ID

NP_078808.3

Blocking

EBP09825

Peptide

Peptide with sequence C-KLHYASYEPTIR, from the internal region of the protein sequence according to

Immunogen NP 078808.3.

Product Comments

This antibody is expected to recognize isoform 1 (NP_078808.3) only.

Peptide Sequence

C-KLHYASYEPTIR

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Shipping Instructions

Refrigerated

Predicted

Human, Mouse, Dog, Cow

Species Reactive

Human **Species**

Human **Gene ID**

79582

Mouse

66722 Gene ID

Product

Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$

IHC Results Paraffin embedded Human Testis. Recommended concentration: 5µg/ml.

ELISA

Antibody detection limit dilution 1:32000. Detection

Limit

Type

Application Pep-ELISA, IHC

GALLERY IMAGES

