

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SPONDIN 2 / MINDIN ANTIBODY

SKU: EB09965

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias SPON2|spondin 2, extracellular matrix protein|M-spondin|Mindin|DIL-1|DIL1

Names

Accession NP 036577.1 ID

Blocking

EBP09965 **Peptide**

Peptide with sequence C-RKNQYVSNGLRDFA, from the internal region of the protein sequence according to Immunogen

NP 036577.1.

Product

Reported variants NP_036577.1 and NP_001121797.1 represent identical protein. Comments

Peptide

C-RKNOYVSNGLRDFA Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human, Mouse, Rat, Cow

Human

10417

Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Limit

Antibody detection limit dilution 1:16000.

Western

Blot

Approx 35kDa band observed in preliminary testing of Human and Mouse Prostate lysate after 1µg/ml antibody staining (calculated MW of 35.8kDa according to Human NP 036577.1). Primary incubation 1 hour at room

temperature.

Application

Pep-ELISA **Type**





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600