

GOAT ANTI-SPONDIN 2 / MINDIN ANTIBODY

SKU: EB09965



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	SPON2 spondin 2, extracellular matrix protein M-spondin Mindin DIL-1 DIL1
Accession ID	NP_036577.1
Blocking Peptide	EBP09965
Immunogen	Peptide with sequence C-RKNQYVSNGLRDFA, from the internal region of the protein sequence according to NP_036577.1.
Product Comments	Reported variants NP_036577.1 and NP_001121797.1 represent identical protein.
Peptide Sequence	C-RKNQYVSNGLRDFA
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Cow
Human Gene ID	10417
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 35kDa band observed in preliminary testing of Human and Mouse Prostate lysate after 1µg/ml antibody staining (calculated MW of 35.8kDa according to Human NP_036577.1). Primary incubation 1 hour at room temperature.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)