

GOAT ANTI-SRC ANTIBODY

SKU: EB05154



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	v-src avian sarcoma (Schmidt-Ruppin A-2) viral oncogene homolog tyrosine-protein kinase SRC-1 tyrosine kinase pp60c-src protooncogene SRC, Rous sarcoma proto-oncogene tyrosine-protein kinase SRC p60-Src c-SRC SRC1 ASV HGNC:11283 v-src sarcoma (Schmidt-Ruppin A-2) viral oncogene homolog (avian) SRC
Accession ID	NP_005408.1; NP_938033
Immunogen	Peptide with sequence GSNKSKPKDASQR-C, from the N Terminus of the protein sequence according to NP_005408.1; NP_938033.
Product Comments	This antibody is expected to recognise both reported isoforms (NP_005408 and NP_93803
Peptide Sequence	GSNKSKPKDASQR-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	6714
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 60kDa band seen in A431 cell lysates. Recommended for use at 0.5-2µg/ml. Overnight incubation is recommended with this antibody.