



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-SRC ANTIBODY**

**SKU:** EB05154

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / v-src avian sarcoma (Schmidt-Ruppin A-2) viral oncogene homolog|tyrosine-protein kinase SRC-1|tyrosine

**Alias** kinase pp60c-src|protooncogene SRC, Rous sarcoma|proto-oncogene tyrosine-protein kinase SRC|p60-Src|c-**Names** SRC|SRC1|ASV|HGNC:11283|v-src sarcoma (Schmidt-Ruppin A-2) viral oncogene homolog (avian)|SRC

**Accession** NP\_005408.1; NP\_938033

Peptide with sequence GSNKSKPKDASQR-C, from the N Terminus of the protein sequence according to **Immunogen** 

NP 005408.1; NP 938033.

GSNKSKPKDASQR-C

**Product** This antibody is expected to recognise both reported isoforms (NP\_005408 and NP\_93803

Comments **Peptide** 

Sequence Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** Refrigerated Instructions

**Predicted** Human, Mouse, Rat, Dog **Species** 

Reactive Human **Species** 

Human 6714 Gene ID

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade

Western Western blot: Approx 60kDa band seen in A431 cell lysates. Recommended for use at 0.5-2µg/ml. Overnight

**Blot** incubation is recommended with this antibody.