

GOAT ANTI-SRCASM / TOM1L1 ANTIBODY

SKU: EB05328



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	OK/KNS-CL.3 target of myb1 (chicken)-like 1 target of myb1 (chicken) homolog-like 1 Src activating and signaling molecule target of myb 1-like 1 SRCASM target of myb1-like 1 (chicken) TOM1L1
Accession ID	NP_005477.2
Blocking Peptide	EBP05328
Immunogen	Peptide with sequence C-IDAQHKGAQNDGD, from the C Terminus of the protein sequence according to NP_005477.2.
Peptide Sequence	C-IDAQHKGAQNDGD
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	10040
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:9000.
Western Blot	No signal obtained yet but low background observed in Human Brain and Kidney lysates at upto 2µg/ml.
Application Type	Pep-ELISA