

GOAT ANTI-SRY ANTIBODY

SKU: EB07265



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	testis-determining factor sex-determining region on Y sex determining region protein essential protein for sex determination in human males TDY TDF sex determining region Y SRY
Accession ID	NP_003131
Blocking Peptide	EBP07265
Immunogen	Peptide with sequence C-HREKYPNYKYRP, from the internal region of the protein sequence according to NP_003131.
Peptide Sequence	C-HREKYPNYKYRP
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	6736
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Western Blot: Approx 38kDa band observed in Human Brain lysates (calculated MW of 23.9kDa according to NP_003131.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Taketo et al, Dev Dyn. 2005 Jun;233(2):612-22.; PMID: 15789443). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB