Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SS-56 / TRIM68 ANTIBODY

SKU: EB05896

SPECIFICATIONS

Unit Size 100 μg

Synonym /

Alias Ro/SSA1 related protein FLJ10369|FLJ10369|SSA protein SS-56|SS-56

Names

Accession NP_060543.5 ID

Blocking

EBP05896 **Peptide**

Immunogen Peptide with sequence TAPLAICSLDGED, from the C Terminus of the protein sequence according to NP_060543.5.

Peptide

Sequence

TAPLAICSLDGED

Shipping

Refrigerated Instructions

Predicted

Human, Cow **Species**

Human

55128 **Gene ID**

Product

Grade

ELISA

Detection

Antibody detection limit dilution 1:4000.

Limit Western

Blot

Preliminary experiments in Human Spleen and Prostate lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

results been reported with other antibodies/lysates?

Application

Pep-ELISA Type