

GOAT ANTI-SS-56 / TRIM68 ANTIBODY

SKU: EB05896



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	Ro/SSA1 related protein FLJ10369 FLJ10369 SSA protein SS-56 SS-56
Accession ID	NP_060543.5
Blocking Peptide	EBP05896
Immunogen	Peptide with sequence TAPLAICSLDGED, from the C Terminus of the protein sequence according to NP_060543.5.
Peptide Sequence	TAPLAICSLDGED
Shipping Instructions	Refrigerated
Predicted Species	Human, Cow
Human Gene ID	55128
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Preliminary experiments in Human Spleen and Prostate lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

