



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SSBP2 ANTIBODY

SKU: EB09712

Telephone: (650) 697-3600





250kDa

150kDa

100kDa 75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Unit Size

100 µg

Synonym / Alias

single-stranded DNA-binding protein 2|HSPC116|DKFZp686F03273|single-stranded DNA binding protein

2|SSBP2 Names

Accession

NP_036578.2

Blocking Peptide

EBP09712

Immunogen

Peptide with sequence YGKGKSNSSAVPSD-C, from the N Terminus of the protein sequence according to

NP 036578.2.

Peptide

Sequence

YGKGKSNSSAVPSD-C

Refrigerated Instructions

Predicted

Species

Human, Mouse, Dog, Cow

Reactive

Species

Human

Human

23635

Gene ID Mouse

66970 **Gene ID**

Rat Gene ID 361877

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade ELISA

Antibody detection limit dilution 1:32000. Detection

Limit

Approx 38kDa band observed in Human Heart lysates (calculated MW of 37.8kDa according to NP_036578.2). Western

Blot Recommended concentration: 0.3-1µg/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa

150kDa

100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa