

GOAT ANTI-SSBP2 ANTIBODY

SKU: EB09712

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	single-stranded DNA-binding protein 2 HSPC116 DKFZp686F03273 single-stranded DNA binding protein 2 SSBP2
Accession ID	NP_036578.2
Blocking Peptide	EBP09712
Immunogen	Peptide with sequence YGKGKSNSAVPSD-C, from the N Terminus of the protein sequence according to NP_036578.2.
Peptide Sequence	YGKGKSNSAVPSD-C
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Dog, Cow
Reactive Species	Human
Human Gene ID	23635
Mouse Gene ID	66970
Rat Gene ID	361877
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 38kDa band observed in Human Heart lysates (calculated MW of 37.8kDa according to NP_036578.2). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa