

GOAT ANTI-STAP2 / BKS (AA391-403) ANTIBODY

SKU: EB05074

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	STAP2 signal transducing adaptor family member 2 BKS FLJ20234 brk kinase substrate BRK substrate breast tumor kinase substrate brk kinase substrate signal-transducing adaptor protein 2 signal-transducing adaptor protein-2
Accession ID	NP_060190.2; NP_001013863.1
Blocking Peptide	EBP05074
Immunogen	Peptide with sequence C-LQKKLEKRRALEH, from the C Terminus of the protein sequence according to NP_060190.2; NP_001013863.1.
Product Comments	This antibody is expected to recognize isoform 1 (NP_060190.2) and isoform 2 (NP_001013863.1).
Peptide Sequence	C-LQKKLEKRRALEH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	55620
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Small Intestine. Recommended concentration: 3.75µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 49kDa band observed in Human Breast lysates (calculated MW of 49.6kDa according to NP_060190.2). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

