


GOAT ANTI-STAT4 (AA 369-383) ANTIBODY

SKU: EB09862



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	STAT4 SLEB11 signal transducer and activator of transcription 4 OTTHUMP00000163553
Accession ID	NP_003142.1
Blocking Peptide	EBP09862
Immunogen	Peptide with sequence TLSNRRFVLCGTNVK, from the internal region of the protein sequence according to NP_003142.1.
Peptide Sequence	TLSNRRFVLCGTNVK
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	6775
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	In transfected HEK293 transiently expressing full-length Human STAT4 (MYC and DYKDDDDK tagged), a band of approx.85kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was observed using anti-myc tag antibody. Recommended concentration, 0.5-1µg/ml.
Application Type	Pep-ELISA, WB-Trf

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

