
Telephone: (650) 697-3600



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

## **GOAT ANTI-STAT4 (AA 369-383) ANTIBODY**

**SKU:** EB09862







Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **SPECIFICATIONS**

**Unit Size** 100 µg

Synonym /

Alias STAT4|SLEB11|signal transducer and activator of transcription 4|OTTHUMP00000163553

Names

**Accession** NP\_003142.1 ID

**Blocking** 

EBP09862 **Peptide** 

Peptide with sequence TLSNRRFVLCGTNVK, from the internal region of the protein sequence according to **Immunogen** 

NP\_003142.1.

**Peptide** 

TLSNRRFVLCGTNVK Sequence

Shipping

Refrigerated Instructions

**Predicted** 

Human **Species** 

Reactive

**Species** 

Human

Human Gene ID

6775

**Product** 

**Grade** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**ELISA** 

**Detection** 

Antibody detection limit dilution 1:1000.

Limit

Western **Blot** 

In transfected HEK293 transiently expressing full-length Human STAT4 (MYC and DYKDDDDK tagged), a band of approx.85kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was

observed using anti-myc tag antibody. Recommended concentration, 0.5-1µg/ml.

**Application** 

Pep-ELISA, WB-Trf **Type** 

## **GALLERY IMAGES**





Telephone: (650) 697-3600



