

GOAT ANTI-STK36 ANTIBODY

SKU: EB09334



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	serine/threonine kinase 36 fused KIAA1278 FU DKFZp434N0223 serine/threonine kinase 36, fused homolog (Drosophila) STK36
Accession ID	NP_056505.1
Blocking Peptide	EBP09334
Immunogen	Peptide with sequence C-QSLPHSSRPASAK, from the C Terminus of the protein sequence according to NP_056505.1.
Peptide Sequence	C-QSLPHSSRPASAK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog
Human Gene ID	27148
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)