

GOAT ANTI-STRA8 (AA186-198) ANTIBODY

SKU: EB11545

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	STRA8 stimulated by retinoic acid gene 8 protein homolog stimulated by retinoic acid gene 8 homolog
Names	(mouse) OTTHUMP00000207988
Accession ID	NP_872295.1
Blocking Peptide	EBP11545
Immunogen	Peptide with sequence C-QNLSEERKASLRQ, from the internal region of the protein sequence according to NP_872295.1.
Peptide Sequence	C-QNLSEERKASLRQ
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	346673
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 35kDa band observed in lysates of cell line HEK293 and Human Testis (calculated MW of 36.9kDa according to NP_872295.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa