

GOAT ANTI-STRA8 (AA70-81) ANTIBODY

SKU: EB11544



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	STRA8 stimulated by retinoic acid gene 8 protein homolog stimulated by retinoic acid gene 8 homolog (mouse) OTTHUMP00000207988
Accession ID	NP_872295.1
Blocking Peptide	EBP11544
Immunogen	Peptide with sequence C-NNLRKTVYSQSD, from the internal region of the protein sequence according to NP_872295.1.
Peptide Sequence	C-NNLRKTVYSQSD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Rat
Human Gene ID	346673
Rat Gene ID	500079
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)