

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-STRA8 (AA70-81) ANTIBODY**

**SKU:** EB11544

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

STRA8|stimulated by retinoic acid gene 8 protein homolog|stimulated by retinoic acid gene 8 homolog **Alias** 

(mouse)|OTTHUMP00000207988 **Names** 

Accession

NP 872295.1 ID

**Blocking** 

EBP11544 **Peptide** 

Peptide with sequence C-NNLRKTVYSQSD, from the internal region of the protein sequence according to Immunogen

NP 872295.1.

**Peptide** 

C-NNLRKTVYSQSD Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

**Species** 

Human, Rat

Human Gene ID

346673

**Rat Gene ID** 500079

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:8000.

Limit

**Application** Pep-ELISA **Type**