


GOAT ANTI-STS2 / TULA (AA470-482) ANTIBODY

SKU: EB10884



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names UBASH3A| ubiquitin associated and SH3 domain containing A| CLIP4| STS-2| TULA| OTTHUMP00000109368| T-cell ubiquitin ligand protein| cbl-interacting protein 4| suppressor of T-cell receptor signaling 2

Accession ID NP_061834.1; NP_001001895.1

Blocking Peptide EBP10884

Immunogen Peptide with sequence C-DSGIRISSVFASP, from the internal region of the protein sequence according to NP_061834.1; NP_001001895.1.

Product Comments This antibody is expected to recognize both reported isoforms (NP_061834.1; NP_001001895.1).

Peptide Sequence C-DSGIRISSVFASP

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	53347
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Breast shows perinuclear staining in the epithelial cells of the tubules. Recommended concentration, 3-5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 70kDa band observed in Human Spleen lysates (calculated MW of 69.8kDa according to NP_001001895.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

