


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-STS2 / TULA (AA470-482) ANTIBODY

SKU: EB10884



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias **Names** UBASH3A| ubiquitin associated and SH3 domain containing A| CLIP4| STS-2| TULA| OTTHUMP00000109368| T-

cell ubiquitin ligand protein| cbl-interacting protein 4| suppressor of T-cell receptor signaling 2

Accession

NP_061834.1; NP_001001895.1

Blocking

ID

EBP10884 **Peptide**

Peptide with sequence C-DSGIRISSVFASP, from the internal region of the protein sequence according to **Immunogen**

NP_061834.1; NP_001001895.1.

Product Comments

This antibody is expected to recognize both reported isoforms (NP_061834.1; NP_001001895.1).

Peptide Sequence

C-DSGIRISSVFASP

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Shipping Instructions Refr

Refrigerated

Predicted

Human

Species Reactive

Human

Species

mannan

Human Gene ID

53347

Product

Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results

In paraffin embedded Human Breast shows perinuclear staining in the epithelial cells of the lubules.

Recommended concentration, 3-5µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:8000.

Limit

Western

Approx 70kDa band observed in Human Spleen lysates (calculated MW of 69.8kDa according to

Blot NP_001001895.1). Recommended concentration: 1-3μg/ml.

Application

Type


Pep-ELISA, WB, IHC

GALLERY IMAGES



Telephone: (650) 697-3600





250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa