



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-STX6 ANTIBODY

SKU: EB07214



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa







SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias HGNC:11441|syntaxin 6|STX6

Names

ID

Accession

NP_005810.1

Blocking

EBP07214 **Peptide**

Peptide with sequence C-QALAERKNRQA, from the internal region of the protein sequence according to **Immunogen**

NP 005810.1.

Peptide

C-QALAERKNRQA Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog **Species**

Reactive

Human **Species**

Human

10228 Gene ID

Mouse

58244 Gene ID

Rat Gene ID 60562

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Approx 30kDa band observed in Human Lymph Node lysates (calculated MW of 29.2kDa according to Western

NP_005810.1). In transfected HEK293 transiently expressing Human STX6 a band of approx. 30kDa is **Blot** observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 0.3-1µg/ml.

Application

Pep-ELISA, WB **Type**

Telephone: (650) 697-3600



Email: customerservice@vectorlabs.com

GALLERY IMAGES

2	e	Ю	lk	D	a
1	5	0	kI	D	а

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

