

GOAT ANTI-STX6 ANTIBODY

SKU: EB07214



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias HGNC:11441|syntaxin 6|STX6

Names

Accession ID NP_005810.1

Blocking Peptide EBP07214

Immunogen Peptide with sequence C-QALAERKNRQA, from the internal region of the protein sequence according to NP_005810.1.

Peptide Sequence C-QALAERKNRQA

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat, Dog

Reactive Species Human

Human Gene ID 10228

Mouse Gene ID 58244

Rat Gene ID 60562

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Limit Antibody detection limit dilution 1:32000.

Western Blot Approx 30kDa band observed in Human Lymph Node lysates (calculated MW of 29.2kDa according to NP_005810.1). In transfected HEK293 transiently expressing Human STX6 a band of approx. 30kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 0.3-1µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

