

GOAT ANTI-SULF1 ANTIBODY

SKU: EB09199

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	sulfatase FP SULF-1 KIAA1077 HSULF-1 FLJ41750 FLJ38022 FLJ30905 sulfatase 1 SULF1
Names	
Accession ID	NP_055985.2
Blocking Peptide	EBP09199
Immunogen	Peptide with sequence C-QSYNKEKGVKKQ, from the internal region of the protein sequence according to NP_055985.2.
Product Comments	Reported variants represent identical protein: NP_055985.2, NP_001121676.1, NP_001121678.1, NP_001121677.1.
Peptide Sequence	C-QSYNKEKGVKKQ
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	23213
Mouse Gene ID	240725
Rat Gene ID	171396
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx. 110kDa band observed in HumanTestis lysates (calculated MW of 101kDa according to NP_055985.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

