


# GOAT ANTI-SWAP70 ANTIBODY

**SKU:** EB05384



## SPECIFICATIONS

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

**Storage Instructions** Aliquot and store at -20°C. Minimize freezing and thawing.

**Synonym / Alias Names** SWAP switching B-cell complex 70kDa subunit|HSPC321|FLJ39540|KIAA0640 protein|KIAA0640|SWAP-70|SWAP-70 protein|SWAP70

**Accession ID** NP\_055870.2; NP\_001284643.1

**Blocking Peptide** EBP05384

**Immunogen** Peptide with sequence C-EEREKNWKEKKTTE, from the C Terminus of the protein sequence according to NP\_055870.2; NP\_001284643.1.

**Peptide Sequence** C-EEREKNWKEKKTTE

**Purification Method** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

**Shipping Instructions** Refrigerated

<b>Predicted Species</b>	Human
<b>Reactive Species</b>	Human
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>
<b>IHC Results</b>	Paraffin embedded Human Tonsil. Recommended concentration: 5µg/ml.
<b>ELISA</b>	
<b>Detection Limit</b>	Antibody detection limit dilution 1:128000.
<b>Application Type</b>	Pep-ELISA, IHC

## DOCUMENTS

- [Data Sheet](#)

## GALLERY IMAGES

