



GOAT ANTI-SYNAPTOGRIN 3 ANTIBODY

SKU: EB06585



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	synaptogyrin 3 MGC:20003 SYNGR3
Usage Summary	Flow Cytometry: A customer reported weak staining on U251 cells at a concentration of 10ug/ml.
Accession ID	NP_004200.2
Blocking Peptide	EBP06585
Immunogen	Peptide with sequence C-QRTAPGPATTQAGD, from the internal region of the protein sequence according to NP_004200.2.
Peptide Sequence	C-QRTAPGPATTQAGD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Rat, Dog
Human Gene ID	9143
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

DOCUMENTS

- [Data Sheet](#)