

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600



GOAT ANTI-SYNAPTOGYRIN 3 ANTIBODY

SKU: EB06585

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias synaptogyrin 3|MGC:20003|SYNGR3

Names

Usage Flow Cytometry: A customer reported weak staining on U251 cells at a concentration of

Summary 10ug/ml.

Accession

NP_004200.2

ID

Blocking EBP06585 Peptide

Peptide with sequence C-QRTAPGPATTQAGD, from the internal region of the protein sequence according to **Immunogen**

NP_004200.2.

Peptide

C-QRTAPGPATTQAGD Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Rat, Dog Species

Human

9143 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Application Pep-ELISA **Type**